11 Publication number:

0 064 680

**A2** 

12)

## **EUROPEAN PATENT APPLICATION**

(21) Application number: 82103677.9

22 Date of filing: 29.04.82

(5) Int. Cl.<sup>3</sup>: C 12 N 15/00 C 12 N 1/20 //C12R1/13, C12R1/15

30 Priority: 30.04.81 JP 65777/81 25.09.81 JP 151464/81

(43) Date of publication of application: 17.11.82 Bulletin 82/46

(84) Designated Contracting States: BE CH DE FR GB IT LI NL SE Applicant: KYOWA HAKKO KOGYO CO., LTD.
 Ohtemachi Bidg., 6-1 Ohtemachi Itchome
 Chiyoda-ku Tokyo 100(JP)

(7) Inventor: Katsumata, Ryoichi Popuragaoka Coopo 6-401, 2-12-3, Naruse Machida-shi Tokyo(JP)

(7) Inventor: Oka, Tetsuo 2360-17, Nara-machi Midori-ku Yokohama-shi Kanagawa-ken(JP)

Inventor: Furuya, Akira 1-10-5, Tamami Tama-ku Kawasaki-shi Kanagawa-ken(JP)

(2) Representative: Vossius Vossius Tauchner Heunemann Rauh Siebertstrasse 4 P.O. Box 86 07 67 D-8000 München 86(DE)

54 Novel lysozyme-sensitive microorganism.

(5) A novel lysozyme-sensitive microorganism belonging to the genus *Corynebacterium* or *Brevibacterium* and having a sensitivity to lysozyme at a concentration of less than 25 µg/ml is provided from selected mutants. This novel microorganism is especially suitable for use in recombinant DNA technology.

April 29, 1982

Kyowa Hakko Kogyo Co., Ltd., Tokyo, Japan "NOVEL LYSOZYME-SENSITIVE MICROORGANISM"

## Background of the Invention

The present invention relates to a novel lysozymesensitive microorganism belonging to the genus <u>Corynebacterium</u> or <u>Brevibacterium</u>, a process for preparing protoplasts
from the microorganism, a process for transforming the microorganism, a process for using the microorganism as a host in
recombinant DNA technology, and the obtained transformants.

Tysozyme-sensitive microorganisms belonging to the genus Corynebacterium or Brevibacterium have been known (Japanese Published Unexamined Patent Application No. 122794/79, etc.). However, the microorganisms of the present invention are more lysozyme-sensitive than those of the Japanese Patent Application, and therefore useful for gene engineering studies as hereinafter described.

20

25

15

5

10

# Detailed Description of the Invention

The present novel lysozyme-sensitive microorganism is a microorganism belonging to the genus <u>Corynebacterium</u> or <u>Brevibacterium</u> and having a sensitivity to lysozyme at a concentration of less than 25 µg/ml. The microorganism of the present invention can be converted to protoplasts treating it only with lysozyme under hypertonic conditions to remove the cell wall, without pretreatment with a drug such as an antibiotic during growth.

30

Recently, the genetic engineering techniques of inserting a foreign gene in plasmid or phage DNA in vitro and introducing the resultant recombinant DNA into host cells have attracted attention.

35

The transformants obtained by this technology are important as a means to produce the foreign DNA by the autonomous replication of the recombinant DNA and to endow the cell with valuable properties due to the presence of the foreign DNA.

Most of the research in this field has been carried out using Escherichia coli as a host, and attempts have also been made to establish "recombinant DNA technology" with other industrially important microorganisms than Escherichia coli, such as Bacillus subtilis capable of producing amylase, etc., actinomyces capable of producing antibiotics, etc., yeasts capable of producing alcohol, etc.

5

10

15

30

In genetic engineering, it is necessary to isolate plasmids for preparation of vectors and analysis of a hybrid plasmid wherein foreign DNA is inserted. To this end cells must be ruptured first. The rupture must be carried out under mild conditions to avoid shearing of DNA. Thus, mechanical rupturing such as ultrasonication and French press is unsuitable, and usually a procedure using a cell wall-lytic enzyme is used. Lytic enzymes suitable for each of the aforesaid microbial species on which research of recombinant DNA is made have been known and are used.

coli or Bacillus subtilis can be readily removed with
lysozyme, i.e. an enzyme capable of lysing cell walls of
bacteria. However, there are some bacteria whose cell walls
are hardly lysed with lysozyme. For these bacteria, such a
procedure has been used that the bacterium is contacted
with a high concentration of an amino acid such as glycine
or with an antibiotic such as penicillin during growth to
render cells lysozyme-sensitive, and the resulting bacterium
is treated with lysozyme.

Corynebacterium or Brevibacterium, which are used in industrial production of various amino acids such as glutamic acid and lysine, usually have cell walls that are resistant to lysozyme. Therefore, such a procedure for removing the cell walls as mentioned above has been so far used for these bacteria species, that is, the cells are treated first with an antibiotic such as penicillin during the growth and then with lysozyme. However, the procedure is complicated and has the additional disadvantage that the regenerating ability of the resulting protoplast is low.

To establish a recombinant DNA technology in useful bacteria such as those belonging to the genus Corynebacterium or Brevibacterium, the present inventors have been investigating processes for removing the cell walls of these bacteria. The present inventors studied a 5 process for removing cell walls with lysozyme, using lysozymesensitive strains of the bacteria species belonging to the genus Corynebacterium or Brevibacterium disclosed in Japanese Published Unexamined Patent Application No. 122794/79, 114293/80 or 42593/81. However, the cell walls of the 10 strains were not removed with lysozyme without pretreatment with penicillin, etc. as in the case of the lysozymeinsensitive strains. As a result of further study, the present inventors have found that the cell walls of lysozymeultrasensitive strains, which have a much higher sensitivity 15 to lysozyme than the lysozyme-sensitive strains disclosed in the above-mentioned Japanese Application, can be lysed with lysozyme without any pretreatment with penicillin, etc. Furthermore, the present inventors have found that plasmids can be easily extracted from these lysozyme-ultrasensitive 20 cells, and the regenerating ability of protoplasts formed by the lysozyme treatment is high, so that the transformation of the protoplasts by DNA can be efficiently carried out. The present invention has been completed on the basis of these findings. 25

The lysozyme-sensitive strains of the present invention can be obtained by subjecting lysozyme-insensitive strains (for example, wild-type strains) or the lysozyme-sensitive strains of the genus Corynebacterium or Brevibacterium to mutation treatment (hereinafter, the lysozyme-sensitive strains of the present invention are referred to as "lysozyme-ultrasensitive strains"). Conventional mutation procedures, for example, ultraviolet irradiation,  $\gamma$ -ray irradiation, chemical mutagenesis and the like can be used. In selecting lysozyme-ultrasensitive strains from mutagenized cells, it is necessary to select those which fail to grow in a medium containing lysozyme at a very low concentration (less than 25  $\mu$ g/ml), in which the parent strains can fully

30

grow, and which can grow in a medium containing no lysozyme. Examples of the present microorganisms include the L-15 strain (FERM P-5946) derived from Corynebacterium glutamicum 225-106 (FERM P-5945), L-103 strain (FERM P-5947) derived from Corynebacterium herculis ATCC 13868, L-204 5 strain (FERM P-5948) derived from Brevibacterium divaricatum ATCC 14020 and L-312 strain (FERM P-5949) derived from Brevibacterium lactofermentum ATCC 13655. These lysozymeultrasensitive strains were deposited with the Fermentation Research Institute and assigned the accession numbers 10 mentioned above, and further deposited with the American Type Culture Collection, USA under the accession numbers ATCC 31834, 31833, 31866, 31867 and 31868, respectively. The 225-106 strain named "Corynebacterium glutamicum"

is a strain isolated from a soil sample and has the following bacteriological properties. Investigation of the bacteriological properties was conducted according to the procedure
described in "Manual of Microbiological Methods" by the
Society of American Bacteriologist Committee on Bacteriological Technique (1957).

- - On an agar plate, a single, circular, lustrous and pale yellow colony; On a slant, a similar pale yellow opaque colony; On an agar stab, abundant growth on surface and slight growth in deep; In a liquid medium, slight growth and slightly flocculent sediment.

# III. Physiological characteristics

30

35

optimum temperature 25 - 37°C; growth occurs slightly at 42°C

		- 5 -	0004880
	2)	рн:	optimum pH 7-8; growth occurs at pH 6-9
	3)	Thermal resistance:	none
	4)	Requirement of free oxygen:	aerobic
5	5)	Gelatin liquefaction:	none
	6)	Assimilation of casein:	negative
•	7)	Indole production:	none
10	8)	Catalase:	positive
	9)	Assimilation of starch:	negative
•			

- 10) Acid production from glucose, fructose, mannose and maltose; non-acid production from xylose, galactose, lactose and glycerol
- 11) Requirement of biotin : positive

15

20

35

- 12) Glutamic acid is accumulated in a large quantity in a medium wherein the amount of biotin is restricted.
- 13) Lactic acid and α-ketoglutaric acid are accumulated in a medium containing biotin in a high concentration.
- 25 These characteristics are compared with those of bacteria disclosed in J. Gen. Appl. Microbiol., 73, 279 301 (1967). Since the characteristics coincide well with those of Corynebacterium glutamicum, the 225-106 strain is identified as a strain of Corynebacterium glutamicum.
- An example of a procedure for obtaining a lysozymeultrasensitive strain of the present invention from a lysozyme-insensitive or sensitive strain of the genus Corynebacterium or Brevibacterium is described below.

A parent strain is inoculated in a nutrient medium and cultured at 20 - 35°C with shaking. Culturing is discontinued in the course of the logarithmic growth phase. The cells are collected, washed with physiological saline

0064680

solution and then suspended in an appropriate buffer solution (pH 5.5-8.0), for example, M/20 Tris-maleate buffer (pH 6.0) in a concentration of 1 x  $10^7$  - 1 x  $10^9$  cells/ml. To the resulting suspension is added nitrosoguanidine in a final concentration of 200-500 µg/ml and the mixture is allowed to stand at 20-35°C for 30 minutes to 1 hour.

5

10

15

20

25

30

35

The cells are collected by centrifugation, washed with the same buffer and suspended in physiological saline solution. A certain amount of the suspension diluted with physiological saline solution is spread on a nutrient agar medium, for example, NB agar medium (medium containing 1.8% agar in NB medium which contains 20 g of bouillon powder and 5 g of yeast extract in 1  $\ell$  of water). After incubation at  $25-35^{\circ}\text{C}$  for 1 to 4 days, the formed colonies are replicaplated to a nutrient agar medium and a nutrient agar medium containing  $2.5-25~\mu\text{g/ml}$  lysozyme. The replica plates are incubated at  $25-35^{\circ}\text{C}$  for 1 to 4 days, and then those colonies which can grow on the nutrient agar medium but not on the nutrient agar medium containing lysozyme are kept as lysozyme-ultrasensitive mutants.

All the aforementioned lysozyme-ultrasensitive strains are obtained in this way.

The lysozyme-sensitivity of the lysozyme-ultrasensitive mutants can be determined by inoculating a drop of cell suspension containing about 10<sup>4</sup> cells of the mutant in the logarithmic phase of growth on NB agar medium containing lysozyme in a series of two-fold dilution and by observing the growth after incubation at 30°C for 2 days.

A minimal lysozyme concentration at which no growth of a given microorganism is observed is taken as a lysozyme sensitivity of the microorganism (minimum inhibitory concentration), and the comparison of sensitivity of the mutants with that of the respective parent strains is given in Table 1.

The lysozyme-ultrasensitive mutants of the present invention are obtained by selecting the strains unable to grow on the NB agar medium containing  $2.5-25~\mu g/ml$  lysozyme,

as described above, whereas the lysozyme-sensitive mutants of bacteria species belonging to the genus Corynebacterium or Brevibacterium disclosed in Japanese Published Unexamined Patent Application No. 122794/79, were obtained by selecting the strains unable to grow on the NB agar medium containing  $200-400 \mu g/ml$  lysozyme. Thus, these two seem to have a possibility of being essentially different in properties. Although the lysozyme sensitivity of the parent strains used in these two experiments as measured by the method described above are almost comparable, between 400 and 800 µg/ml, that of the resultant mutants differsconsiderably. That is, the sensitivities of the mutants of the present invention is 1.6-3.2 µg/ml, whereas that of the mutants of the Japanese Patent Application is 25 - 400 µg/ml. Furthermore, as will be explained below, the lysozyme-ultrasensitive mutants of 15 the present invention have advantages due to the high lysozyme susceptibility of their cell walls, which can be lysed without pretreatment with penicillin, whereas the lysozyme-sensitive mutants of the Japanese Patent Application have no such characteristics. Thus, it is apparent that the lysozymeultrasensitive strains are distinct from the lysozymesensitive strains.

25

20

5

10

30

#### Table 1

5	<u>Strains</u>	FERM P	ATCC.	vity, minimum inhi- bitory concentra- tion (MIC, µg/ml)
	Corynebacterium glutamicum 225-106	5945	31833	800
10	Corynebacterium glutamicum L-15	5946	31834	3.2
	Corynebacterium herculis ATCC 13868	· •		400
	Corynebacterium herculis L-103	5947	31866	3.2
15	Brevibacterium divaricatum ATCC 14020	· <b>-</b>	٠	800
	Brevibacterium divaricatum L-204	5948	31867	1.6
20	Brevibacterium lactofermentum ATCC 13655	-		400
	Brevibacterium lactofermentum L-312	5949	31868	3.2

The usefulness of the present microorganisms
whose cell walls can be easily lysed and removed by lysozyme
resides in that 1) efficient transformation of protoplasts
by DNA is possible due to the high regenerating ability of
the protoplasts and 2) extraction of DNA such as plasmid DNA
from the cells can be easily made.

The transformation and the extraction are described in detail below.

1) Preparation of protoplasts from the lysozymeultrasensitive mutants and transformation of the protoplasts by DNA

5

10

15

20

25

30

35

Protoplasts can be prepared by treating the cultured cells only with lysozyme. Any medium wherein the inventive microorganisms propagate may be used. For example, NB medium can be used. Microorganisms are inoculated onto a medium and cultured at 25-37°C with shaking up to the middle to

late in the logarithmic phase of growth. The cells are collected and suspended in an appropriate hypertonic medium containing 0.2-10 mg/ml lysozyme. An example of the hypertonic medium is PFM medium containing 0.4 M sucrose and 0.01 M MgCl<sub>2</sub>·6H<sub>2</sub>O in two-fold diluted SSM medium and being adjusted to pH 7.0 - 8.5. The SSM medium consists of 20 g of glucose, 10 g of  $(NH_4)_2SO_4$ , 3 g of urea, 1 g of yeast extract, 1 g of  $KH_2PO_4$ , 0.4 g of  $MgCl_2 \cdot 6H_2O$ , 10 mg of  $FeSO_4 \cdot 7H_2O$ , 0.2 mg of MnSO<sub>4</sub>·4-6H<sub>2</sub>O, 0.9 mg of ZnSO<sub>4</sub>·7H<sub>2</sub>O, 0.4 mg of CuSO<sub>4</sub>·5H<sub>2</sub>O, 0.09 mg of Na<sub>2</sub>B<sub>4</sub>O<sub>7</sub>·10H<sub>2</sub>O, 0.04 mg of (NH<sub>4</sub>)<sub>6</sub>Mo<sub>7</sub>O<sub>24</sub>·4H<sub>2</sub>O, 30 µg of biotin, and 1 mg of thiamine hydrochloride in 1 £ of water and is adjusted to pH 7.2. Normal cells are gradually converted to protoplasts during incubation at 30 - 37°C, and the degree of the conversion can be monitored under an optical microscope. The time necessary for the conversion of most of the cells to protoplasts depends upon the strain used, and usually is 0.5 - 8 hours under the above-mentioned conditions. After the lysozyme treatment, the number of normal cells resistant to hypotonic condition is less than 10<sup>-4</sup> of the initial number of the normal cells subjected to the lysozyme treatment.

The protoplasts thus prepared have the ability to form colonies (regenerating ability) on a hypertonic agar medium. An example of a hypertonic agar medium is a medium containing 13 g/2 agar in RCGP medium [a medium containing

0064680

5 g of glucose, 5 g of casein hydrolyzate, 2.5 g of yeast extract, 3.5 g of  $K_2HPO_4$ , 1.5 g of  $KH_2PO_4$ , 0.41 g of MgCl<sub>2</sub>.  $6H_{2}O$ , 10 mg of  $FeSO_{4} \cdot 7H_{2}O$ , 2 mg of  $MnSO_{4} \cdot 4 - 6H_{2}O$ , 0.9 mg of  $z_{n}s_{0_{4}}\cdot 7_{H_{2}}$ 0, 0.04 mg of  $(NH_{4})_{6}Mo_{7}O_{24}\cdot 4H_{2}O$ , 30 µg of biotin, 2 mg of thiamine hydrochloride, 135 g of disodium succinate, and 30 g of polyvinylpyrrolidone (molecular weight: 10,000) in 1 % of water and being adjusted to pH 7.2]. of protoplastSon .RCGP agar medium depends upon the strain used, and the concentration and time of lysozyme treatment, and the efficiency of regeneration is usually 20 - 70% on the basis of the initial number normal cells subjected of to the lysozyme treatment. The time required for observable regeneration of colonies depends upon the strain used, and it usually takes 5 - 7 days for the colonies to become large enough to pick up.

Transformation of the protoplasts obtained from the microorganisms having the above-mentioned characteristics is carried out by mixing the protoplasts with DNA under hypertonic conditions and in the presence of polyethyleneglycol (PEG) or polyvinyl alcohol (PVA), and divalent metal cations.

Transformants can be obtained as normal cells by regenerating the protoplasts in or on a hypertonic agar medium such as RCGP agar medium in the same manner as in the regeneration of the protoplasts mentioned above. Transformants can be selected on the basis of phenotypic properties which genes on the introduced DNA provide. The selection can be made in or on a hypertonic agar medium simultaneously with the regeneration, or after the regeneration.

In the transformation of the protoplasts by plasmid 30 DNAs according to the foregoing procedure, transformants can be obtained with a frequency of  $10^{-4} - 10^{-6}$  per regenerated cell.

Transformation is described in detail in the examples shown below.

5

10

15

20

## 2) Isolation of plasmid DNA from transformant

The microorganism is cultured in a medium suitable for propagation of the microorganism, for example, NB medium, collected and treated with lysozyme to lyse the cell walls. Plasmidscan be easily isolated from the lysate according to the ordinary procedure as described, for example, in Biochim. Biophys. Acta 383 457 - 463 (1975).

Suitable concentration of lysozyme depends upon the strain used, and is usually 0.1-10 mg/ml.

5

30 -

35

Application of recombinant DNA technique to the 10 bacteria species of the genera Corynebacterium and Brevibacterium, which are used in the industrial production of. various amino acids such as glutamic acid, lysine, etc., has such a possibility that genes taking part in biosynthesis of useful substances such as amino acid, etc. from these 15 bacteria and other microorganisms can be cloned to intensify the biosynthetic activities by the amplification of the genetic information, thereby increasing the production of useful substances. Also the genes of animals and plants can be cloned, so that useful polypeptides can be produced in these 20 bacteria species by the expression of the genetic informa-Thus the application has tremend ous industrial poten-The lysozyme-ultrasensitive microorganisms of the tial. present invention enable extractive separation of DNA such as plasmid, etc. from the cells and transformation of the 25 cells by DNA, and therefore can facilitate application of the recombinant DNA technique to the genera Corynebacterium and Brevibacterium.

The present inventors have further found that the microorganisms of the present invention are useful as safe. host microorganisms in the recombinant DNA technology.

It is prerequisite for recombinant DNA techniques to establish a safe host-vector system. This requires very low pathogenicity, no parasitic relationship to humans, and rapid loss of viability when released into the environment.

Escherichia coli, the first host/in this technology, is by nature a parasite to humans, and its safety has long been

discussed. However, at least strain K-12 cultured in laboratories for a long period of time was found to have considerably reduced parasitic tendencies and pathogenicity and has been approved as B l host. It is reported that, when the K-12 strain was administered to humans, the microorganism could be recovered from the body over a period of a few days with an overall recovery of about 1% [Science 209 391-394 (1980)].

In extending the recombinant DNA technique to microorganisms capable of producing amino acids, nucleic acids, antibiotics, etc., it is necessary to establish safe host-vector system. Most of the relevant microorganisms are isolated from soils and thus their parasitic relationship to humans is regarded as low, but some may survive the conditions prevailing in the digestive tract of the human body, though not colonizing it, and might be recovered alive. To increase the safety of such strains, the property of perishing in the digestive tract can be artificially endowed on the strains. In the case of a B 2 host, such a weak property as perishing when released into the environment has been endowed.

It has been found that the lysozyme-ultrasensitive microorganisms of the present invention perish much more readily in the digestive tracts of animals, when tested therein. That is, the microorganisms shown in Table 2 were orally administered at a dose of about 3-5 x 10 per male ddy mouse (body weight: 19 g ± 1), each group consisting of two mice. The number of microorganisms remaining in the intestines and the number of microorganisms discharged with feces were measured by spreading appropriately diluted 10% homogenates each thereof on NB plate medium. As is apparent from Table 2, the lysozyme-ultrasensitive microorganisms of the present invention perish in a considerable degree, and are discharged, as compared with the parent strains, and the number of the microorganisms remaining in the intestines is also more rapidly decreased.

Table 2

	Microorganisms used	Recovery %					
5		In feces		In intestines			
		* hours 0 - 20.5	hours 20.5-43.5	** hours 20.5	hours 43.5		
	Corynebacterium glutamicum 225-106	53.8	0.08	5.7 x 10 <sup>-7</sup>	<9 x 10 <sup>-10</sup>		
10	Corynebacterium glutamicum . L-15	4.3	<9.3 x 10 <sup>-9</sup>	<1.3 x 10 <sup>-9</sup>	$< 8.7 \times 10^{-10}$		
	Corynebacterium herculis ATCC 13868	44.5	0.05	$4.6 \times 10^{-7}$	<8.5 x 10 <sup>-10</sup>		
15 .	Corynebacterium herculis L-103	3.4	$< 8.7 \times 10^{-9}$	<1 x 10 <sup>-9</sup>	$<7.3 \times 10^{-10}$		
	Corynebacterium divaricatum ATCC 14020	57.2	0.09	6.4 x 10 <sup>-7</sup>	<8.7 x 10 <sup>-10</sup>		
20	Corynebacterium divaricatum L-204	7.4	<8.9 x 10 <sup>-9</sup>	$<1.1 \times 10^{-9}$	<9.1 x 10 <sup>-10</sup>		
	Brevibacterium lactofermentum ATCC 13655	64.5	0.1	$4.3 \times 10^{-7}$	<8.2 x 10 <sup>-9</sup>		
2.5	Brevibacterium lactofermentum L-312	5.6	$<7.4 \times 10^{-6}$	$<8.4 \times 10^{-9}$	<2.3 x 10 <sup>-9</sup>		

\* Time after the administration of the relevant microorganism. Column "0-20.5" shows the recovery of
microorganisms (% of the given microorganisms) in
total feces from zero hour (the time of administration)
up to 20.5 hours.

30

35

\*\* Column "20.5 hours" shows the recovery of microorganisms in the intestine removed from mice at 20.5 hours. The present inventors have studied the sensitivities of a lysozyme-ultrasensitive strain toward various drugs. A drop of the cell suspension corresponding to about 10<sup>4</sup> cells in the logarithmic phase of growth is inoculated on an NB agar medium containing a test drug in a two-fold serial

dilution. After incubation at 30°C for 2 days, the growth of each inoculum was scored. As shown in Table 3, it was found that the strain L-15 is more sensitive to penicillin, rifampicilin and of course to lysozyme than its parent 225-106. This suggests that these lysozyme-ultrasensitive strains are more readily killed or inhibited from growing than their parents, and therefore the antibiotic resistance conferred to host cells by plasmid vectors will be counterbalanced.

15

10

5

Table 3

20	

	Sensitivity to drugs (*MIC, µg/ml)		
	lysozyme	penicillin G	rifampicilin
Corynebacterium glutamicum 225-106	800	0.12	0.016
Corynebacterium glutamicum L-15	3.2	0.015	0.004

25

30

As described above, the present lysozyme-ultrasensitive microorganisms are considerably easier to destroy in the digestive tract of animals, and are highly sensitive to various drugs. This shows that the present microorganisms are preferable as host microorganisms for recombinant DNA.

### Example 1

35

Corynebacterium glutamicum 225-106, Corynebacterium herculis ATCC 13868, Brevibacterium divaricatum ATCC 14020, and Brevibacterium lactofermentum ATCC 13655 are inoculated in NB medium (pH 7.2), and cultured at 30°C with shaking.

<sup>\*</sup> MIC: Minimum inhibitory concentration

The growth is terminated in the middle of the logarithmic phase of growth, and cells are collected, washed with physiological saline solution, and then suspended in M/20 Trismaleate buffer (pH 6.0) to make about  $5 \times 10^8$  cell/ml.

The suspensions are mixed with nitrosoguanidine to make a final concentration of 400 µg/ml, and allowed to stand at 25°C for 30 minutes. The cells are collected by centrifugation, washed with the same buffer, and then suspended in physiological saline solution. The suspension is appropriately diluted with the saline and spread on NB-agar medium. Incubation is carried out at 30°C for 2 days and the resultant colonies are replica-plated to NB-agar medium with or without 12.5 µg/ml lysozyme (egg white, 6 x crystallized, made by Seikagaku Kogyo Co., Ltd.). The replica-plates are incubated at 30°C for 2 days and then those colonies which can grow on the NB medium but not on the NB medium containing lysozyme are kept as lysozyme-ultra-

20

sensitive mutants.

15

5

10

Table 4

Parent strains	Lysozyme-ultrasensitive mutant strains		
Corynebacterium glutamicum 225-106	Corynebacterium glutamicum L-15		
Corynebacterium herculis ATCC 13868	Corynebacterium herculis L-103		
Brevibacterium divaricatum ATCC 14020	Brevibacterium divaricatum L-204		
Brevibacterium lactofermentum ATCC 13655	Brevibacterium L-312		
	Corynebacterium glutamicum 225-106  Corynebacterium herculis ATCC 13868  Brevibacterium divaricatum ATCC 14020  Brevibacterium		

#### 35 Example 2

Lysis and removal of cell walls of the microorganisms of the present invention by lysozyme and regeneration of the protoplasts formed.

Seed cultures of Corynebacterium glutamicum 225-106, Corynebacterium herculis ATCC 13868, Brevibacterium devaricatum ATCC 14020, Brevibacterium lactofermentum ATCC 13655, and the lysozyme-ultrasensitive mutant strains obtained in Example 1 by derivation therefrom are respectively inoculated onto NB medium and cultured at 30°C with shaking. Absorbancy (OD) at 660 nm is measured with a colorimeter made by Tokyo Koden Co. A small portion of the culture is appropriately diluted with physiological saline solution when OD reaches 0.6. An aliquot of the diluted solution is spread onto NB agar medium to determine the cell number. At the same time cells are harvested from the culture broth and suspended in a medium (pH 7.6) containing 1 mg/ml lysozyme in RCGP medium to make the cell density about 109 The suspension is transferred to an L-shaped test tube and incubated at 30°C with gentle shaking.

Five hours thereafter, cells are recovered by centrifugation at 2,500 x G for 10 minutes, suspended in PFM medium, twice washed by centrifugation and resuspended in PFM medium to prepare protoplast suspension. A portion of the suspension is appropriately diluted with RCGP medium, and spread onto RCGP agar medium, while another portion is appropriately diluted with physiological saline solution, and spread onto NB agar medium. They are cultured at 30°C to determine the number of colony-forming cells under hypotonic condihypertonic conditions as well as under tions The number of colonies formed under hypotonic conditions is measured on the second day of culture (the increase with further incubation). number does not number of colonies formed under the hypertonic conditionsis measured on the seventh day of the culture after which increases no further. the number of regenerated colonies The results are shown in Table 5.

30

10

15

20

Table 5

	Strains	Number of lyso- zyme-treated	Number of colony-forming cells after lysozyme treatment		
5		cells (cfu/ml)	Hypotonic	Hypertonic	
•	Corynebacterium glutamicum L-15	1.8 x 10 <sup>9</sup>	9.8 x 10 <sup>3</sup>	6.7 x 10 <sup>8</sup>	
10	Corynebacterium glutamicum 225-106	2.3 x 10 <sup>9</sup>	2.1 x 10 <sup>9</sup>	2.3 x 10 <sup>9</sup>	
	Corynebacterium herculis L-103	1.2 x 10 <sup>9</sup>	6.5 x 10 <sup>4</sup>	2.4 x 10 <sup>8</sup>	
15	Corynebacterium. herculis ATCC 13868	1.3 × 10 <sup>9</sup>	1.3 x 10 <sup>9</sup>	1.4 x 10 <sup>9</sup>	
	Brevibacterium divaricatum L-204	9.6 x 10 <sup>8</sup>	1.1 x 10 <sup>3</sup>	6.4 x 10 <sup>8</sup>	
20	Brevibacterium divaricatum ATCC 14020	1.2 × 10 <sup>9</sup>	1.3 x 10 <sup>9</sup>	1.4 x 10 <sup>9</sup>	
	Brevibacterium lactofermentum L-312	1.0 x 10 <sup>9</sup>	1.4 × 10 <sup>4</sup>	5.3 x 10 <sup>8</sup>	
25	Brevibacterium lactorermentum ATCC 13655	1.1 x 10 <sup>9</sup>	9.9 x 10 <sup>8</sup>	1.0 x 10 <sup>9</sup>	
4.5		<u>                                       </u>	1	<u> </u>	

With all the parent strains, no formation of spherical protoplasts is observed and only normal cells of the short rod form are observed under a microscope after they are subjected to action of lysozyme. The number of colonies formed under hypotonic conditions is not substantially different before and after the lysozyme treatment. On the other hand, the lysozyme-ultrasensitive strains are converted to protoplasts by the action of lysozyme, and the ratio of the cells capable of forming colonies under hypotonic conditions where the protoplast is killed by rupture (osmotic-resistant normal cells) is

30

10<sup>-4</sup> - 10<sup>-6</sup> per cell subjected to the lysozyme treatment. These protoplasts can form colonies under hypertonic conditions, and their regeneration efficiency is 20 - 70% of the number of the cells subjected to lysozyme treatment. Cells from the regenerated colonies are normal, so they can grow on NB agar medium.

#### Example 3

5

10

15

20

25

30

35

Transformation of the protoplasts prepared from the microorganisms of the present invention by plasmid pCG4 DNA

## Preparation of plasmid pCG4

Plasmid pCG4 for use in transformation in the present Example is extracted from the cultured cells of Corynebacterium glutamicum 225-250 (FERM P-5939), a microorganism carrying pCG4, in the following manner.

Five milliliters of seed culture of Corynebacterium glutamicum 225-250 in NB medium is inoculated into 400 ml of SSM medium. Culturing is carried out at 30°C with shaking. Optical density (OD) at 660 nm is measured with a Tokyo Koden colorimeter and at OD 0.2 penicillin G is added to the broth in a final concentration of 0.5 U/ml. Further, growth is continued at 30°C to OD about 0.6.

Cells are recovered from the culture broth, washed with TES buffer (pH 8.0) consisting of 0.03M tris(hydroxymethyl)aminomethane (Tris), 0.005M EDTA and 0.05M NaCl and suspended in a lysozyme solution (pH 8.0) consisting of 25% sucrose, 0.1M NaCl, 0.05M Tris and 0.8 mg/ml lysozyme to make 10 ml of a suspension. The suspension is incubated at Then, 2.4 ml of 5M NaCl, 0.6 ml of 0.5M 37°C for 4 hours. EDTA (pH 8.5) and 4.4 ml of a solution containing 4 g sodium laurylsulfate and 0.7M NaCl are added successively to the reaction suspension. After the mixture is stirred slowly, it is kept on an ice-water bath for 15 hours. The whole lysate is put into a centrifugation tube and centrifuged under 69,400 x g at 4°C for 60 minutes to obtain a super-To the supernatant is added 10% by weight natant.

of polyethyleneglycol 6,000. After the mixture is stirred slowly to dissolve it, it is kept on an ice-water bath. After 16 hours, the mixture is subjected to centrifugation under 1,500 x g for 10 minutes to obtain a pellet. pellet is redissolved in 5 ml of TES buffer and 2.0 ml of 1.5 mg/ml ethidium bromide is added. Cesium chloride is added to the mixture to adjust the density to 1.580. solution is subjected to centrifugation under 105,000 x q at 18°C for 48 hours. After the density gradient centrifugation, a circular DNA closed with covalent bond is detected by UV irradiation as a high density band located in the lower part of the centrifugation tube. The band is taken out from the side of the tube with an injector to obtain a fraction containing/plasmid pCG4. To remove ethidium bromide, the fraction is treated five times with an equal amount of cesium chloride saturated isopropyl alcohol solution consisting of 90% by volume of isopropyl alcohol and . -10% TES buffer. Then, the residue is subjected to dialysis against TES buffer. Thus, 15 µg of plasmid pCG4 is obtained. Plasmid pCG4 is a novel plasmid having a molecular

Plasmid pCG4 is a novel plasmid having a molecular weight of about 19 megadaltons, carrying genes for resistance to streptomycin and spectinomycin and having 4, 7, 9, 6 and 6 restriction sites for Eco RI, Bam HI, Hind III, Pst I and Sal I respectively.

25

30

35

20

5

10

The lysozyme treated cells prepared in the same manner as in Example 2 are transformed with pCG4. Zero point five milliliters of protoplast suspension is taken into a small test tube, centrifuged at 2,500 x G for 5 minutes. The precipitate is resuspended in 1 ml of TSMC buffer (10 mM magnesium chloride, 30 mM calcium chloride, 50 mM Tris, 400 mM sucrose: pH 7.5) and washed by centrifugation. The precipitated protoplasts are resuspended in 0.1 ml of TSMC buffer with gentle shaking.

To the suspension is added 0.1 ml of a DNA solution containing 0.2 µg of DNA which is prepared by diluting a pCG4 solution two-fold. with two-fold concentrated TSMC

buffer and the resultant mixture is stirred. Then, 0.8 ml of a solution containing 20% polyethylene glycol 6,000 in TSMC buffer is added to the mixture and the resultant mixture is gently stirred. Three minutes thereafter, the suspension is mixed with 2 ml of RCGP medium, and centrifuged at 2,500 x G for 5 minutes to remove the supernatant. The precipitated protoplasts are suspended in 1 ml of RCGP medium, diluted to several levels with RCGP medium, and spread onto an RCGP agar medium containing 400 µg/ml of spectinomycin. To determine the number of colony-forming cells, the protoplast suspension is highly diluted, spread onto an RCGP agar. medium without spectinomycin and incubated at 30°C for 6 days.

The colonies formed on the RCGP medium containing spectinomycin are plated onto an NB agar medium containing 12.5 µg/ml of streptomycin and cultured at 30°C for 2 days to investigate cross resistance to spectinomycin and streptomycin. The strains showing cross resistance are selected at random, and/plasmid is isolated according to the procedure shown in Example 4.

15

20

The results are summarized in Table 6.

Table 6

5	Recipient	Frequency of appearance of	Properties of spectino- mycin-resistant cells **		
	microorganism	spectinomycin- resistant cells *	Streptomycin cross resistance	Presence of pCG4	
	Corynebacterium glutamicum 225-106	1.7 x 10 <sup>-7</sup>	- (5/5)	- (2/2)	
10	Corynebacterium glutamicum L-15	$2.4 \times 10^{-3}$	÷ (50 / 50)	+ (2/2)	
	Corynebacterium herculis ATCC 13868	<1.3 x 10 <sup>-7***</sup>			
15	Corynebacterium herculis L-103	$7.5 \times 10^{-4}$	+ (50 / 50)	+ (2/2)	
	Brevibacterium divaricatum ATCC 14020	3.5 x 10 <sup>-7</sup>	- (5/5)	- (2/2)	
20	Brevibacterium divaricatum L-204	1.0 x 10 <sup>-6</sup>	+ (50 / 50)	+ (2/2)	
25	Brevibacterium lactofermentum ATCC 13655	$3.1 \times 10^{-7}$	- (5/5)	- (2/2)	
	Brevibacterium lactofermentum L-312	2.4 x 10 <sup>-6</sup>	+ (50 / 50)	+ (2/2)	

\* Frequency per colony-forming microorganism (regenerated microorganism) under hypertonic conditions.

30

35

\*\*\* Spectinomycin-resistant cells are not obtained from the strain.

The spectinomycin-resistant strains obtained by treating lysozyme-insensitive strains with lysozyme and pCG4

<sup>\*\*</sup> The number of cells tested (denominator) and those of negative (-) or positive (+) cells (numerator) are illustrated in parentheses.

plasmid show no streptomycin cross resistance, and pCG4 is not detected in these spectinomycin-resistant strains by the isolation procedure used for the strain 225-250. On the other hand, the spectinomycin-resistant strains obtained by treating lysozyme-ultrasensitive strains with lysozyme and pCG4 plasmid have/streptomycin cross resistance, and the plasmid isolated from these strains give the same DNA fragments upon digestion with restriction endonuclease Hind III as those from pCG4.

The FERM P numbers of the transformants thus obtained are shown in Table 7.

pCG4 15 Recipient strain PERM P ATCC transformant L15 / pCG4 5950 Corynebacterium 31835 glutamicum L15 L103 / pCG4 Corynebacterium 5951 31869 <u>L10</u>3 herculis 20 Brevibacterium L204 / pCG4 31870 5952 divaricatum L204 L312 / pCG4 Brevibacterium 5953 31871 lactofermentum L312

Table 7

25

30

35

5

10

#### Example 4

Isolation of pCG4 plasmid from the cells transformed with pCG4

Five milliliters of seed culture in NB medium of each pCG4 transformant obtained in Example 3, Corynebacterium glutamicum L15/pCG4, Corynebacterium herculis L103/pCG4,

Brevibacterium divaricatum L204/pCG4 or Brevibacterium lactofermentum L312/pCG4 is inoculated in 400 ml of NB medium and cultured with shaking at 30°C for 12 hours.

The cells are collected from the culture broth, washed with TES buffer, suspended in 10 ml of a lysozyme

solution (25% sucrose, 0.1M NaCl, 0.05M Tris, 0.8 mg/ml lysozyme: pH 8) and subjected to incubation at 37°C for one hour. Then the reaction solution is successively admixed with 2.4 ml of 5M NaCl, 0.6 ml of 0.5M EDTA (pH 8.5) and 4.4 ml of a solution containing 4% sodium laurylsulfate and 0.7M NaCl, gently mixed, and placed in ice-water for 15 hours. The whole lysate is transferred into a centrifugation tube, and subjected to centrifugation at 69,400 x G at 4°C for 60 minutes to obtain a supernatant fraction.

10 Polyethylene glycol 6000 is added to the supernatant to make the final concentration 10% (W/W) with gentle stirring and the resultant solution is placed in ice water. Sixteen hours thereafter, pellets are recovered by centrifugation at 1,500 x G for 10 minutes, and dissolved in 5 ml of TES buffer.

15 The solution is admixed with 2 ml of 1.5 mg/ml ethidium bromide and then with cesium chloride to give a density of 1.580. The solution is subjected to density gradient centrifugation under the same conditions as described in Example 3 to separate plasmid. The fraction containing pCG4 is treated with the isopropyl alcohol solution in the same manner as in Example 3 to extract and remove ethidium bromide, and dialyzed against the TES buffer.

In this manner,  $15-30~\mu g$  of pCG4 DNA was obtained from all the lysozyme-ultrasensitive strains containing pCG4.

# What is claimed is:

5

10

15

20

25

30

- (1) A novel lysozyme-sensitive microorganism belonging to the genus <u>Corynebacterium</u> or <u>Brevibacterium</u> and having a sensitivity to lysozyme at a concentration of less than 25 µg/ml.
- (2) A microorganism according to claim 1, wherein the microorganism is selected from Corynebacterium glutamicum, Corynebacterium herculis, Brevibacterium divaricatum and Brevibacterium lactofermentum.
- (3) A microorganism according to claim 2, wherein the microorganism is selected from Corynebacterium glutamicum L-15 (FERM P-5946, ATCC 31834), Corynebacterium herculis L-103 (FERM P-5947, ATCC 31866), Brevibacterium divaricatum L-204 (FERM P-5948, ATCC 31867) and Brevibacterium lactofermentum L-312 (FERM P-5949, ATCC 31868).
- organism belonging to the genus <u>Corynebacterium</u> or <u>Brevibacterium</u> which comprises treating a microorganism belonging to the genus <u>Corynebacterium</u> or <u>Brevibacterium</u> and having a sensitivity to lysozyme at a concentration of less than 25 µg/ml with lysozyme to form a protoplast of the microorganism.
- (5) A process for transforming a microorganism belonging to the genus <u>Corynebacterium</u> or <u>Brevibacterium</u> which comprises treating a microorganism belonging to the genus <u>Corynebacterium</u> or <u>Brevibacterium</u> and having a sensitivity to lysozyme at a concentration of less than 25 µg/ml with lysozyme, introducing deoxyribonucleic acid into the resulting protoplast and regenerating the protoplast.
- (6) A process for using a microorganism belonging to the genus Corynebacterium or Brevibacterium and having a sensitivity to lysozyme at a concentration of less than 25 µg/ml as a host in recombinant DNA technology.

- (7) A novel lysozyme-sensitive microorganism belonging to the genus <u>Corynebacterium</u> or <u>Brevibacterium</u>, having a sensitivity to lysozyme at a concentration of less than 25 µg/ml and carrying pCG4 plasmid.
- (8) A microorganism according to claim 7, wherein the microorganism is selected from Corynebacterium glutamicum, Corynebacterium herculis, Brevibacterium divaricatum and Brevibacterium lactofermentum.

5

(9) A microorganism according to claim 8, wherein the microorganism is selected from Corynebacterium glutamicum

L15/pCG4 (FERM P-5950, ATCC 31835), Corynebacterium herculis
L103/pCG4 (FERM P-5951, ATCC 31869), Brevibacterium

divaricatum L204/pCG4 (FERM P-5952, ATCC 31870) and Brevibacterium lactofermentum L312/pCG4 (FERM P-5953, ATCC 31871).